

Technical Data Sheet

GENERAL INFORMATION		
PRODUCT NAME REFERENCE	Thermolabile dsDNase I expressed in yeast, lyophilized MT01U-L1DNATXA	
PRODUCTS PROVIDED	01U-L1DNATXA Thermolabile ds	DNase I, lyophilized
UNITS EXPRESSION SYSTEM PURITY	1,000 U Komagataella phaffii ≥ 90%	
DESCRIPTION	Thermolabile double strand DNase I is an endonuclease that cleaves phosphodiester bonds in DNA to release oligonucleotides with 5'-phosphorylated and 3'-hydroxylated ends. Thermolabile dsDNase I has a particularly strong preference for double-stranded DNA. In the presence of Mg ²⁺ as only divalent cation and using oligos as substrate, the activity towards ssDNA is minimum compared to dsDNA. This is why the enzyme can be used to specifically degrade dsDNA, leaving ssDNA and RNA intact. This enzyme is easily inactivated by heat treatment at 55 °C and it can be irreversible inactivated adding DTT at a final concentration of 1 mM while heating. Lyophilized dsDNase I is a freez-dried version of its well- characterized liquid equivalent.	
DELIVERY CONDITION		
DELIVERY CONDITION 01U-L1DNATXA	1,000 U of freeze-dried thermolab	ile dsDNase I
	1,000 U of freeze-dried thermolab	ile dsDNase I
01U-L1DNATXA	1,000 U of freeze-dried thermolab A) For long-term: 20 mM Tris glycerol. B) For short-term: 20 mM Tris	pH 7.5, 10 mM MgCl ₂ , 50 %
01U-L1DNATXA RELEVANT INFORMATION SUGGESTED	A) For long-term: 20 mM Tris glycerol.	рН 7.5, 10 mM MgCl₂, 50 % pH 7.5, 10 mM MgCl₂
01U-L1DNATXA RELEVANT INFORMATION SUGGESTED RECONSTITUTION BUFFERS SUGGESTED 10X REACTION	A) For long-term: 20 mM Tris glycerol.B) For short-term: 20 mM Tris	pH 7.5, 10 mM MgCl₂, 50 % pH 7.5, 10 mM MgCl₂ ²
01U-L1DNATXA RELEVANT INFORMATION SUGGESTED RECONSTITUTION BUFFERS SUGGESTED 10X REACTION BUFFER	 A) For long-term: 20 mM Tris glycerol. B) For short-term: 20 mM Tris 200 mM Tris pH 7.5, 100 mM MgCl 1- Mix the reaction mixture or COMPONENTS 	pH 7.5, 10 mM MgCl₂, 50 % pH 7.5, 10 mM MgCl₂ ²
01U-L1DNATXA RELEVANT INFORMATION SUGGESTED RECONSTITUTION BUFFERS SUGGESTED 10X REACTION BUFFER	 A) For long-term: 20 mM Tris glycerol. B) For short-term: 20 mM Tris 200 mM Tris pH 7.5, 100 mM MgCl 1- Mix the reaction mixture or COMPONENTS RNA* 	pH 7.5, 10 mM MgCl ₂ , 50 % pH 7.5, 10 mM MgCl ₂ ² n ice: 50 μL REACTION ~ 1 μg
01U-L1DNATXA RELEVANT INFORMATION SUGGESTED RECONSTITUTION BUFFERS SUGGESTED 10X REACTION BUFFER	 A) For long-term: 20 mM Tris glycerol. B) For short-term: 20 mM Tris 200 mM Tris pH 7.5, 100 mM MgCl 1- Mix the reaction mixture or COMPONENTS RNA* Sample 	pH 7.5, 10 mM MgCl ₂ , 50 % pH 7.5, 10 mM MgCl ₂ ² n ice: 50 μL REACTION ~ 1 μg Up to 44 μL
01U-L1DNATXA RELEVANT INFORMATION SUGGESTED RECONSTITUTION BUFFERS SUGGESTED 10X REACTION BUFFER	 A) For long-term: 20 mM Tris glycerol. B) For short-term: 20 mM Tris 200 mM Tris pH 7.5, 100 mM MgCl 1- Mix the reaction mixture or COMPONENTS RNA* Sample dsDNase I Reaction Buffer (10X) 	pH 7.5, 10 mM MgCl ₂ , 50 % pH 7.5, 10 mM MgCl ₂ ² n ice: 50 μL REACTION ~ 1 μg Up to 44 μL 5 μL
01U-L1DNATXA RELEVANT INFORMATION SUGGESTED RECONSTITUTION BUFFERS SUGGESTED 10X REACTION BUFFER	 A) For long-term: 20 mM Trisglycerol. B) For short-term: 20 mM Tris200 mM TrispH 7.5, 100 mM MgCl 1- Mix the reaction mixture or COMPONENTS RNA* Sample dsDNase I Reaction Buffer (10X) dsDNase I 	pH 7.5, 10 mM MgCl ₂ , 50 % pH 7.5, 10 mM MgCl ₂ ² n ice: 50 μL REACTION ~ 1 μg Up to 44 μL 5 μL 1 μL (2 U)
01U-L1DNATXA RELEVANT INFORMATION SUGGESTED RECONSTITUTION BUFFERS SUGGESTED 10X REACTION BUFFER	 A) For long-term: 20 mM Tris glycerol. B) For short-term: 20 mM Tris 200 mM Tris pH 7.5, 100 mM MgCl 1- Mix the reaction mixture or COMPONENTS RNA* Sample dsDNase I Reaction Buffer (10X) 	pH 7.5, 10 mM MgCl ₂ , 50 % pH 7.5, 10 mM MgCl ₂ ² n ice: 50 μL REACTION ~ 1 μg Up to 44 μL 5 μL 1 μL (2 U) To 50 μL

2- Incubate at 37 °C for 30 minutes.

3- Heat inactivate at 55 °C for 15 minutes.



Optional: add DTT at a final concentration of 1 mM for an irreversible inactivation.

ACTIVITY UNIT DEFINITION One Unit of activity is defined as an increase in absorbance at 260 nm of 0.001 per minute at 25 °C on the assay conditions (50 μ g/mL calf thymus DNA in buffer 20 mM Tris, 10 mM MgCl₂, pH 7.5).

STORAGEThe lyophilized product can be handled and store at Room
Temperature at least 12 months. Once reconstituted, it is
recommended to store the solution from -20 °C to -80 °C.
Storage at 4 °C is possible for short term. Avoid multiple
freeze/thaw cycles by storing multiple aliquots at -80 °C.

HEALTH AND SAFETYConsult the Safety Data Sheet for information regarding
hazards and safe handling practises.

QUALITY CONTROL

DNase ACTIVITY ASSAY

Once reconstituted, thermolabile dsDNase I activity is measured for each lot by incubating dsDNase I with calf thymus DNA. For that purpose, 40 U of dsDNase I is incubated with 0.05 mg/mL of thymus DNA at 25 °C and the release of dsDNA I is monitored at 260 nm. The resulting units are then compared with the theoretical units, with an accepted 8% deviation from reference units.

TECHNICAL SUPPORT

If you have any questions, feel free to contact us at <u>hello@levprot.com</u>

Polígono Industrial La Charluca, Calle C, Parcela 10-11 50300 Calatayud, Zaragoza (Spain) Tel. (+34) 976198404 <u>hello@levprot.com</u> www.levprot.com

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

DATE 05/09/2023 VERSION 00